



CAPSULE SIZES	000	00	0	1	2	3	4	5
<b>Weight</b>								
Milligrams	163	118	96	76	61	48	38	28
Tolerance (mg)	±10	±7	±6	±5	±4	±3	±3	±2
<b>Capacity</b>								
Capsule Volume (mL)	1.37	0.91	0.68	0.50	0.37	0.30	0.21	0.13
<b>Powder Density - Capsule Capacity (mg)</b>								
0.6 g/mL	822	546	408	300	222	180	126	78
0.8 g/mL	1096	728	544	400	296	240	168	104
1 g/mL	1370	910	680	500	370	300	210	130
1.2 g/mL	1644	1092	816	600	444	360	252	156
<b>Length of the Capsule Parts (Body and Cap)</b>								
Body (in)	0.874	0.796	0.726	0.654	0.601	0.535	0.480	0.366
Tolerance (in)	±0.018	±0.018	±0.018	±0.018	±0.018	±0.018	±0.018	±0.016
Body (mm)	22.20	20.22	18.44	16.61	15.27	13.59	12.19	9.30
Tolerance (mm)	±0.46	±0.46	±0.46	±0.46	±0.46	±0.46	±0.46	±0.40
Cap (in)	0.510	0.462	0.422	0.385	0.352	0.318	0.284	0.244
Tolerance (in)	±0.018	±0.018	±0.018	±0.018	±0.018	±0.018	±0.018	±0.016
Cap (mm)	12.95	11.74	10.72	9.78	8.94	8.08	7.21	6.20
Tolerance (mm)	±0.46	±0.46	±0.46	±0.46	±0.46	±0.46	±0.46	±0.40
<b>External Diameter*</b>								
Body (in)	0.376	0.322	0.289	0.261	0.239	0.219	0.199	0.184
Body (mm)	9.55	8.18	7.34	6.63	6.07	5.57	5.05	4.68
Cap (in)	0.390	0.336	0.300	0.272	0.250	0.229	0.209	0.193
Cap (mm)	9.91	8.53	7.64	6.91	6.35	5.82	5.32	4.91
<b>Overall Closed Length</b>								
Inches	1.029	0.917	0.854	0.765	0.709	0.626	0.563	0.437
Tolerance (in)	±0.012	±0.012	±0.012	±0.012	±0.012	±0.012	±0.012	±0.016
Millimeters	26.1	23.3	21.7	19.4	18.0	15.9	14.3	11.1
Tolerance (mm)	±0.3	±0.3	±0.3	±0.3	±0.3	±0.3	±0.3	±0.4

\* All tolerances ±0.002 in or ±0.06 mm