



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5700



Head Cover, Bouffant Cap, Circular 24", Medisca Safe-Sense™

DESCRIPTION: Bouffant caps (head covers) are manufactured from strong, lightweight, and breathable spun-bounded polypropylene (SPP) fabric, which protects against dirt, grime, and certain dry particulates in controlled environments. These circular bouffant caps come packaged open, as seen in the product image.

FEATURES:

- ✔ Polypropylene nonwoven material
- ✔ Circular: Single elastic for user comfort and stability with high breathability
- ✔ Elastic head band

TECHNICAL SPECIFICATIONS:

Material Composition	PP
Weight (gsm)	12
Stretch (cm)	61
Color	White

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM