



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 5771**



**Gloves, Blue Nitrile Powder-Free, Low Derma, L - 9"- 5 mil, Non-Sterile, Safe-Sense™**

**DESCRIPTION:** Produced using **LOW DERMA™** technology, these single-use revolutionary new-generation nitrile gloves provide exceptional comfort and superior barrier protection for non-sterile hazardous compounding, administration, disposal and other related chemotherapy activities. They are resistant to puncture, abrasion and snag. These gloves are also very commonly used for non-hazardous compounding.

**FEATURES:**

- ✔ Improves user and product protection
- ✔ Meets ASTM D6319 and tested in accordance to ASTM D5151 – freedom from holes – AQL1.5
- ✔ Meets ASTM D6978: Resistant to permeation by chemotherapy drugs\*
- ✔ Meets FDA regulations for food handling applications
- ✔ Meets ISO 11193-1, the recognized standard by Health Canada
- ✔ Beaded cuff to reduce glove roll down and ripping
- ✔ Chlorinated glove for easy donning and removal
- ✔ Produced with soft polymer nitrile formulation
- ✔ Textured fingertips to provide grip and sensitivity
- ✔ Manufactured in an ISO 13485 certified facility

\*Chemical permeation and chemotherapy drug resistance data will be provided upon request. Users are advised to review the material safety data sheet for then chemicals being used to determine the required level of protection.

**TECHNICAL SPECIFICATIONS:**

Grade	Medical Examination
Material	Nitrile
Surface	Smooth with textured fingers
Color	Blue

Size	Length (mm)	Palm Width (mm)	Cuff Thickness (mm)	Finger Thickness (mm)	Palm Thickness (mm)	Tensile Strength (MPa)		Elongation (%)	
						Before Aging	After Aging	Before Aging	After Aging
S	Min 240	85 ± 5	Min 0.06	0.13 ± 0.03	Min 0.08	14	14	500	400
M		95 ± 5				14	14	500	400
L		106 ± 5				14	14	500	400
XL		116 ± 5				14	14	500	400