



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5773



Isolation Gown w/ Sleeves, Yellow, Large, Non-Sterile, Safe-Sense™

DESCRIPTION: These non-surgical isolation gown provides protection against chemicals, liquids and infective agents.

FEATURES:

- ✔ PPE category III
- ✔ AAMI level 2
- ✔ Full back closure, belt attachment and knitted cuffs
- ✔ Isolation gown with sleeves, neck ties and waist ties
- ✔ Laminated non-woven polyethylene
- ✔ Protective against chemicals – EN 14325:2004
- ✔ Protective against infective agents – EN 14126:2003+AC:2004
- ✔ Liquid barriers performance – ANSI/AAMI PB70:2012
- ✔ Latex-free
- ✔ Low-linting
- ✔ Manufactured in an ISO 13485 facility

*Bouffant cap, goggle, face mask and shoe covers are not included.

TECHNICAL SPECIFICATIONS:

| | |
|---------------|---|
| Fabric | Polypropylene and Polyethylene film 42 g/m ² |
| Knitted Cuffs | 100% Polyester |
| Color | Yellow |