



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5789



Balance, Ohaus Ranger 3000 Bench Scale, 30000 x 10g NTEP, w/Ext. Cal.

DESCRIPTION: The multipurpose Ranger 3000 has all of the capabilities to produce precise results rapidly for a variety of applications. With nine application modes, Ranger 3000 is the ideal scale to meet many industrial weighing needs. Portable and equipped with a rechargeable battery, this scale is beneficial in every corner of the factory and is equipped with functions and features that elevate the Ranger 3000 to a class of its own.

FEATURES:

- ✔ The big, bright display makes Ranger 3000 extremely easy to use
- ✔ Produces accurate results within one second
- ✔ Rapid stabilization virtually eliminates your wait time and allows you to move at a pace determined by you and not your weighing equipment
- ✔ The included easy-to-use & install weigh-below hook offers the functionality to perform specific gravity tests
- ✔ With its convenient integrated carrying handle and rechargeable lead acid battery that can provide power for more than 210 hours between recharges

TECHNICAL SPECIFICATIONS:

Maximum Capacity (lb/kg)	60/30
Readability (g)	1.0
Platform Size, W x L (mm/in)	300 x 225/11.8 x 8.9
Battery Life	210 hours with rechargeable battery
Net Weight (lb/kg)	9.2/4.2
Dimensions, L x W x H (mm/in)	311 x 327 x 117/12.2 x 12.9 x 4.6
Working Environment	14°F – 104°F, 85%RH, non-condensing
Display	LCD with white LED backlight + 3 LED
Platform Construction	Stainless Steel
Linearity (g)	± 2
Repeatability (g)	2

WARRANTY: 2 years

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM