



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5790



**Balance, Ohaus PR Series, 820 x 0.01g (LFT off), NTEP, w/Ext. Cal.
PR822N/E AM**

DESCRIPTION: This PR balance can be connected to point-of sale systems or printers via its standard RS-232 port, or via USB or Ethernet through an interface kit. Featuring an easy-to read, brightly-lit display and a simple user interface, the PR is well-constructed for long-term use, and delivers competitive performance.

FEATURES:

- ✔ Equipped with three essential weighing modes, RS-232 connectivity for data transfer and storage, and internal calibration, making it ideal for routine weighing applications.
- ✔ Equipped with an easy-to-read, bright backlit display and a simple user interface.
- ✔ Easy to operate, with almost no training required.
- ✔ The PR's small footprint saves desktop space while providing a large weighing surface.
- ✔ The PR is durably constructed, and features a stainless steel pan to withstand day-to-day use in the workplace.

TECHNICAL SPECIFICATIONS:

| | |
|------------------------------|------------------------------------|
| Maximum Capacity (lb) | 820 |
| Readability (g) | 0.01 |
| Repeatability (g) | 0.001 |
| Linearity (g) | 0.02 |
| Stabilization Time (sec.) | 1 |
| Pan Diameter Size (in/mm) | 7.1 / 180 |
| Temperature Range | 50°F to 86°F/10°C to 30°C |
| Net Weight (lbs/kg) | 3.5 / 7.7 |
| Dimensions, L x W x H(in/mm) | 12.5 x 7.9 x 3.7 / 3.17 x 201 x 93 |

WARRANTY: 2 years

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM