



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5896



pH Accessories, Electrode (Ohaus ST350)

DESCRIPTION: Use in combination with Starter bench and portable meters. OHAUS Starter Series Electrode: ST350 3-in-1 pH Electrodes have the temperature probe in them. The diameter of the electrode should be smaller than the tube.

FEATURES:

- 3-in-1 refillable pH electrode, with integrated temperature probe
- Suitable for normal samples
- Constructed of glass shaft and built to withstand daily use
- Include a temperature sensor, which powers automatic temperature compensation and ensures accurate measurements
- Fit perfectly in the electrode holders on OHAUS bench meters and electrode clips on all OHAUS portable meters

TECHNICAL SPECIFICATIONS:

Measurement Range	0.00-14 pH
Internal Reference Material	Ag/AgCl
Refilling Reference Electrolyte	3M KCl solution
Working Environment	50°F – 104°F/ 10°C – 40°C, 85%RH, non-condensing
Net Weight (g)	70.0
Temperature Range	0°C – 100°C
Shaft Dimension	120 x 12/4.7 x 0.47
Shaft Material	Glass
Cable Length (m)	1
Temperature Sensor	30 kΩ

WARRANTY: 2 years

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM