



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5921



Data Logger Kit, Binder (0°C – 350°C)

DESCRIPTION: The Binder Data Logger kits for temperature can record temperature data of Binder equipment.

FEATURES:

- ✔ Includes calibration certificate
- ✔ Temperature sensor PT 100 with 2m (6.6 ft) extension cable, plus cable bushings and a sensor mounting bracket for the Binder unit
- ✔ Mounting bracket for attaching the Data Logger onto the Binder unit

TECHNICAL SPECIFICATIONS:

Memory Capacity	64000 values
Interval of measurement	1s to 3h
Measuring range	0°C – 350°C / 32°F–662°F
Dimensions, W x H x D (mm/in)	110 x 85 x 35/4.3 x 3.4 x 1.4
Power Supply	1 x lithium battery 3.6V

APPLICABLE PRODUCTS:

- ✔ 6727: Oven, Forced-Air Convection, Binder APT.line, 4.1
- ✔ 6758: Oven, Forced-Air Convection, Binder Avantgarde.line, 2.1 cu.ft
- ✔ 7840: Convection Oven, Mechanical High Temperature, 1.9 cu.ft.
- ✔ 5920: Oven, Forced-Air, Binder Avantgarde.line 26.2 cu.ft

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM