



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 5931**



## pH Pen Meter, 0.1pH Resolution (Ohaus)

**DESCRIPTION:** Pen meters offer a convenient way to measure pH, or the relative acidity or alkalinity of a sample at a given temperature. The pH pen meters can be completely operated with just one hand, freeing up your other hand to handle sample. Pen meters are the economical option when you are simply looking for meters that will provide accurate measurement without restriction.

### FEATURES:

- ✓ Liquid crystal display (LCD) display
- ✓ ABS Housing
- ✓ Easily replaceable electrodes, automatic shut off feature

### TECHNICAL SPECIFICATIONS:

Measurement Range	0.00 – 14 pH
Accuracy	± 0.1 pH
Measurement Resolution	0.1 pH
Display	LCD
Automatic Shutdown	Yes, after 6 minutes of non-use
Battery	4 AG13 1.5V (micro alkaline batteries)
Battery Life	200 hours with disposable batteries
Dimensions, H x L x W (in/mm)	7.28 x 1.5 x 1.77/185 x 38 x 45
Net Weight (g/oz)	110/3.88
Protection Rating	IP67
Working Environment	50°F – 104°F, 85%RH, non-condensing (10°C – 40°C, 85%RH, non-condensing)
Wrist Strap	Yes

### ADDITIONAL ACCESSORIES:

- ✓ 7070: pH Accessories, pH/ORP Storage Solution, 230mL **(REQUIRED)**

**WARRANTY:** 1 year

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)