



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5964



Balance, Ohaus Adventurer, 120g x 0.0001g, w/AutoCal

DESCRIPTION: Ohaus Adventurer AX124 delivers stability, accuracy and fast weighing results. Weigh samples up to 120 grams with sensitivity to 0.0001 gram graduations. This balance incorporates all of the applications necessary for routine weighing and measurement activities. With a color touchscreen, GLP/GMP compliance capabilities, two USB ports, and ingenious draftshield.

FEATURES:

- ✔ Fast stabilization and reliable operation are enhanced by the AutoCal™ option to ensure solid weighing performance and accurate, repeatable results for standard laboratory applications.
- ✔ The wide viewing angle color touchscreen provides access to all of Adventurer's applications, while dual USB ports with GLP/GMP tracking capabilities assist in the monitoring and reporting of data.
- ✔ The two-piece, top-mounted draftshield doors reduce the footprint of the doors in the rear of the balance when open, and open wide on both sides in order to provide ample room to place samples
- ✔ Two USB ports (one on the front and one on the back) provide connectivity with both a PC and a portable external device at the same time
- ✔ Illuminated level indicator located on the front panel makes adjustment easy in any lighting condition
- ✔ Advanced weighing functions include: part counting, percent, dynamic, density, check weighing, totalization, and formulation modes

TECHNICAL SPECIFICATIONS:

| | |
|-------------------------------|--------------------------------------|
| Dimensions, W x L x H (mm/in) | 230 x 354 x 340/9.05 x 13.93 x 13.39 |
| Maximum Capacity (g) | 120 |
| Calibration | Internal |
| Readability (mg) | 0.1 |
| Repeatability (mg) | 0.1 |

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

| | |
|------------------------|-------------------------------------|
| Pan Size (mm/in) | 90/3.5 |
| Pan Construction | Stainless steel |
| Output | RS-232/USB |
| RS-232 Interface | Yes |
| Stabilization time (s) | 3 |
| Net weight (kg/lbs) | 5.1/11.3 |
| Display | Full-Color 4.3 in Touchscreen WQVGA |

WARRANTY: 2 years

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.