



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5971



Head Cover, Bouffant Cap, Pleated, 21", Safe-Sense™

DESCRIPTION: Bouffant caps (head covers) are manufactured from strong, lightweight, and breathable spun-bounded polypropylene (SPP) fabric, which protects against dirt, grime, and certain dry particulates in controlled environments. These bouffant caps come packaged pleated.

FEATURES:

- ✔ Provides full hair coverage improving worker and product protection
- ✔ 12 gms Polypropylene
- ✔ Pleated: Wide elastic for comfort fit and improved stability
- ✔ Lightweight and breathable material
- ✔ Not made with natural rubber latex
- ✔ Comfortable and secure fit
- ✔ Protection from dry particulates
- ✔ Fluid resistant material for optimal protection
- ✔ Low-linting
- ✔ Disposable
- ✔ White
- ✔ Manufactured in an ISO 13485 facility
- ✔ Meets FDA requirements for food handling applications
- ✔ Non-sterile