



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5973



Gloves, Black, Nitrile Powder-Free, S - 9"- 5 mil, Non-Sterile, Safe-Sense™

DESCRIPTION: These single use gloves provide maximum protection for non-sterile hazardous compounding, administration, disposal and other related chemotherapy activities. They are resistant to puncture, abrasion and snag. These gloves are also very commonly used for non-hazardous compounding.

FEATURES:

- ✔ Improves user and product protection
- ✔ Robust design for superior durability
- ✔ Polyurethane coating for easy donning and removal
- ✔ Beaded cuff to reduce glove roll down and ripping
- ✔ Textured fingertips for exceptional grip and sensitivity
- ✔ Produced with distilled water for a silky feeling
- ✔ Not made with natural rubber latex
- ✔ Anti-static properties
- ✔ Low levels of extractables and particles
- ✔ Meets and exceeds ASTM D5151, ASTM D6978 and ASTM D6319
- ✔ Meets FDA regulations for Food Handling applications
- ✔ Manufactured in an ISO 13485 certified facility

TECHNICAL SPECIFICATIONS:

Grade	Medical Examination
Material	Nitrile
Surface	Smooth with textured fingers
Color	Black

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

Size	Length (mm)	Palm Width (mm)	Cuff Thickness (mm)	Finger Thickness (mm)	Palm Thickness (mm)	Tensile Strength (MPa)		Elongation	
						Before Aging	After Aging	Before Aging	After Aging
S	240 ± 10	83 ± 5	0.07 ± 0.03	0.11 ± 0.03	0.08 ± 0.03	14	14	500	400
M		95 ± 5				14	14	500	400
L		104 ± 5				14	14	500	400
XL		112 ± 5				14	14	500	400

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.