



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 5976



Gloves, Black, Nitrile Powder-Free, XL - 9"- 5 mil, Non-Sterile, Safe-Sense™

DESCRIPTION: These single use gloves provide maximum protection for non-sterile hazardous compounding, administration, disposal and other related chemotherapy activities. They are resistant to puncture, abrasion and snag. These gloves are also very commonly used for non-hazardous compounding.

FEATURES:

- ✔ Improves user and product protection
- ✔ Robust design for superior durability
- ✔ Polyurethane coating for easy donning and removal
- ✔ Beaded cuff to reduce glove roll down and ripping
- ✔ Textured fingertips for exceptional grip and sensitivity
- ✔ Produced with distilled water for a silky feeling
- ✔ Not made with natural rubber latex
- ✔ Anti-static properties
- ✔ Low levels of extractables and particles
- ✔ Meets and exceeds ASTM D5151, ASTM D6978 and ASTM D6319
- ✔ Meets FDA regulations for Food Handling applications
- ✔ Manufactured in an ISO 13485 certified facility

TECHNICAL SPECIFICATIONS:

Grade	Medical Examination
Material	Nitrile
Surface	Smooth with textured fingers
Color	Black

Size	Length (mm)	Palm Width (mm)	Cuff Thickness (mm)	Finger Thickness (mm)	Palm Thickness (mm)	Tensile Strength (MPa)		Elongation	
						Before Aging	After Aging	Before Aging	After Aging
S	240 ± 10	83 ± 5	0.07 ± 0.03	0.11 ± 0.03	0.08 ± 0.03	14	14	500	400
M		95 ± 5				14	14	500	400
L		104 ± 5				14	14	500	400
XL		112 ± 5				14	14	500	400