



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 5999**



## Beaker, Glass w/Spout, 15000 mL (Duran)

**DESCRIPTION:** This reusable glass beaker with heavy wall offers excellent mechanical strength and durability, while providing reduced thermal shock resistance, ideal for use when performing titration.

### FEATURES:

- Made with Borosilicate glass
- Large labelling field for easy marking in fired-on, highly durable, white ceramic
- Spout for clean pouring

### TECHNICAL SPECIFICATIONS:

Height (mm)	390
Capacity (mL)	15,000
Type	Beaker
Material	Borosilicate Glass

USA | CANADA | AUSTRALIA PACIFIC | EUROPE  
[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)