



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6000



Beaker, Glass w/Spout, 300mL

DESCRIPTION: This reusable glass berzelius beaker features a convenient tall-form design. It offers excellent mechanical strength and durability, while providing high resistance to chemical attack and thermal shock, ideal for use when performing titration. It meets ASTM E-960 Type III specifications for classification, design, capacity & dimensions, and markings.

FEATURES:

- ☑ Slightly flared top, reinforcing bead and pouring spout
- ☑ Durable matte finish marking area
- ☑ White scales to approximate volumes
- ☑ Made with Borosilicate glass

TECHNICAL SPECIFICATIONS:

Diameter (in/cm)	2.51/6.4
Height (in/cm)	4.64/11.8
Capacity (mL/oz)	300/10.14
Graduated	Yes
Graduation range (mL)	25 to 250
Increments (mL)	25
Scale	Dual
Disposable	No
Type	Beaker
Style	Berzelius
Material	Borosilicate Glass

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM