



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6007



Applicator, UnoDose, 42ml, 0.25ml per Click, Grey

DESCRIPTION: UnoDose™ Metered Topical Applicator redefines metered-dose applicator technology for optimum provider efficiency and patient ease of use.

FEATURES:

- ✔ In-container electronic mortar and pestle compounding, no transfer required
- ✔ Multiple exit holes for reliable priming
- ✔ Fast prime feature quickly removes air with no repetitive motion, saving time and labor
- ✔ Open, screw-free chamber for easy filling
- ✔ Allows cream application without hand contact
- ✔ Convex, five hole cap for uniform topical application
- ✔ Audible, tactile, and visual metered dose indicators for accurate dosing
- ✔ Locking cover secures against accidental spills
- ✔ Clearly visible reorder indicator
- ✔ Registered as an FDA Class I Medical Device
- ✔ Medicament contacting materials: polyolefins free of phthalates and bisphenol A (BPA)
- ✔ UV Protection: UV-inhibiting materials used for the medicament container

TECHNICAL SPECIFICATIONS:

Dose size	25mL per click (quarter turn of dosing knob)
Maximum usable volume (mL)	42
Residual volume (mL)	1.5

INCLUDES:

PK/25: 4 mixing paddles, 1 mixing lid and 1 priming tool

PK/100: no mixing paddles, no mixing lids and no prime tools (sold separately)

PK/250: no mixing paddles, no mixing lids and no prime tools (sold separately)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM

Please note that the information contained herein is for informational purposes only. The information is based on recognized compendia, manufacturers' instructions and/or other authoritative sources and texts. MEDISCA takes no responsibility for the validity or accuracy of this information. Customer MUST refer to USP, FDA, Health Canada, TGA and/or their State Board to ensure they comply with the specific regulations and ensure that appropriate procedures are in place.