




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6008



Applicator, UnoDose, 42ml, 0.25ml per Click, Purple

DESCRIPTION: UnoDose™ Metered Topical Applicator redefines metered-dose applicator technology for optimum provider efficiency and patient ease of use.

FEATURES:

- ✔ In-container electronic mortar and pestle compounding, no transfer required
- ✔ Multiple exit holes for reliable priming
- ✔ Fast prime feature quickly removes air with no repetitive motion, saving time and labor
- ✔ Open, screw-free chamber for easy filling
- ✔ Allows cream application without hand contact
- ✔ Convex, five hole cap for uniform topical application
- ✔ Audible, tactile, and visual metered dose indicators for accurate dosing
- ✔ Locking cover secures against accidental spills
- ✔ Clearly visible reorder indicator
- ✔ Registered as an FDA Class I Medical Device
- ✔ Medicament contacting materials: polyolefins free of phthalates and bisphenol A (BPA)
- ✔ UV Protection: UV-inhibiting materials used for the medicament container

TECHNICAL SPECIFICATIONS:

Dose size	25mL per click (quarter turn of dosing knob)
Maximum usable volume (mL)	42
Residual volume (mL)	1.5

INCLUDES:

- ✔ PK/25: 4 mixing paddles, 1 mixing lid and 1 priming tool
- ✔ PK/100: 16 mixing paddles, 4 mixing lids and 4 prime tools
- ✔ PK/250: 40 mixing paddles, 10 mixing lids and 10 prime tools

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM