



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6020**



## **Shoe Cover, Non-Skid, Non-Sterile, Safe-Sense™**

**DESCRIPTION:** These disposable anti-skid shoe covers protect shoes from dirt and grime while tread on the bottom provides traction. The spunbond polypropylene construction reduces fluid penetration.

### **FEATURES:**

- ✔ Size: 40cm x 16cm
- ✔ Color: Blue
- ✔ Single elastic band and rim provides non-binding comfort fit and stability
- ✔ Thread pattern reduce risk of slipping
- ✔ Elastic opening for a snug fit
- ✔ Sonic bonded for durability
- ✔ Lightweight and breathable
- ✔ Protection from dry particulates
- ✔ Fluid resistant
- ✔ Not made with natural rubber latex
- ✔ Low-linting
- ✔ Disposable
- ✔ Manufactured in an ISO 13485 facility
- ✔ Meets FDA requirements for food handling applications
- ✔ Non-sterile