



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6033



Weigh Boat, Anti-Static, 3.5 x 3.5 x 1", Medium, White

DESCRIPTION: Anti-static disposable polystyrene weight boats are flexible for handling liquids and solids. Square with stable rounded bottom, smooth surface and contoured sides.

FEATURES:

- ✓ Easy pour spout
- ✓ Boats are shallow with a wide, flat bottom to resist tipping
- ✓ Rounded corners simplify transfer

TECHNICAL SPECIFICATIONS:

| | |
|-------------------|-------------------------------|
| Capacity (mL) | 100 |
| Dimension (mm/in) | 89 x 89 x 25/3.5 x 3.5 x 1 |
| Color | White |
| Material | Polystyrene (PS), Anti-Static |