



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6035**



## **Weigh Canoe, Anti-Static, 2.1 x 1.5 x 0.5", Small, White**

**DESCRIPTION:** This weigh canoe protects the pans of weighing balances from damage and provides a convenient way to transfer weighed materials to another container.

### **FEATURES:**

- ✔ Pour-spout design easily dispenses samples without flexing the dish
- ✔ Flat-bottom dishes are from 12-mil thick polystyrene for better strength, stability and handling

### **TECHNICAL SPECIFICATIONS:**

Capacity (mL)	15
Dimension (mm/in)	53 x 38 x 13/2.1 x 1.5 x 0.5
Color	White
Material	Polystyrene (PS), Anti-Static