



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6045**



## **Weigh Canoe, Anti-Static, 2.7 x 2 x 0.6", Medium, Black**

**DESCRIPTION:** This weigh canoe protects the pans of weighing balances from damage and provides a convenient way to transfer weighed materials to another container.

### **FEATURES:**

- ☑ Pour-spout design easily dispenses samples without flexing the dish
- ☑ Flat-bottom dishes are from 12-mil thick polystyrene for better strength, stability and handling

### **TECHNICAL SPECIFICATIONS:**

Capacity (mL)	25
Dimension (mm/in)	68 x 50 x 16/2.7 x 2 x 0.6
Color	Black
Material	Polystyrene (PS), Anti-Static