



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6106






SPARE KIT PART, U-HEAT

For TTS-8 Tube Sealer (0.25")

DESCRIPTION: The dual element spare parts kit is used for replacement of heating elements and PTFE after prolonged use.

INCLUDED PARTS:

-  2 x heating element
-  2 x PTFE insulation tape 8"
-  4 x 6 mil 2" wide x 10.5" long PFE Fluoropolymer

APPLICABLE PRODUCTS:

-  7050-01: Tube Sealer, 8" Table Top w/Built-In Compressor

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM