



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6119



Stopper, 20 mm, 2 Prongs, Chlorobutyl, Grey Sterile

DESCRIPTION: This all-purpose design for lyophilization applications reduces the possibility of legs sticking to each other. The stopper prongs lock into the interior of the vial neck, providing a tight seal and keeping the stopper from backing out.

FEATURES:

- ✔ Manufactured from the highest quality raw materials
- ✔ Components manufactured to yield exceptional dimensional stability
- ✔ Free of latex, nitrosamines and 2-MCBT (2-mercaptobenzothiazole)
- ✔ Very well suited for ready-for-sterilization and ready-for-use applications
- ✔ Autoclavable

TECHNICAL SPECIFICATIONS:

- ✔ Durometer: 46
- ✔ Style: 2-prong lyophilisation
- ✔ Rubber material: Grey Chlorobutyl

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM