



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6132



Image is for illustrative purposes only. Actual product may differ in color, size, and appearance.

PIONEER SERIES BALANCE, OHAUS (220 × 0.0001 g)

Model PX224, Internal Calibration

DESCRIPTION: The newly designed Pioneer Series balance is used for various essential functions, in a pharmaceutical laboratory, such as basic weighing, parts counting, percent weighing, animal weighing and density determination.

Included in Package: Balance, Power Adapter and Plug, Stainless Steel Pan, Warranty card.

FEATURES:

- ✔ Static removal bar provides a convenient grounding location
- ✔ Large viewing angle with bright backlit LCD display
- ✔ 2 line display provides additional information or guidance prompts to the user
- ✔ InCal mechanism which calibrates the balance automatically without the use of an external calibration mass

WARRANTY: 2 years

TECHNICAL SPECIFICATIONS:

| | |
|--|--|
| Capacity (g) | 220 |
| Readability (mg) | 0.1 |
| Dimensions, W × D × H (cm / in) | 20.9 × 32.1 × 30.9 / 8.2 × 12.6 × 12.2 |
| Platform size, Ø (cm / in) | 9.0 / 3.5 |
| Draft Shield Dimensions, W × D × H (cm / in) | 17.2 × 13.1 × 20.5 / 6.8 × 5.2 × 8.1 |
| Net Weight (Kg / lb) | 4.5 / 9.9 |
| Shipping Weight (Kg / lb) | 7.0 / 15.4 |
| Calibration | Internal, Semi-Automatic |
| Communication | RS232, USB |
| Power Supply | Input: 100 – 240 V, ~200 mA 50 – 60 Hz, 12 – 18 VA Output: 12 VDC, 0.5 A |
| Compliance | IEC / EN 61010-1, CAN / CSA C22.2 61010; UL 61010-1 IEC / EN 61326-1 Class B, Basic Environments; FCC Part 15 Class A; Canada ICES-003 Class A |