



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6146



Stopper, 13 mm, Snap-On, Chlorobutyl, Grey (Sterile)

DESCRIPTION: This design provides maximum sealing for vials with a blow-back glass finish, which is a small ridge on the interior of the vial's neck. The snap-on stopper locks into the blow-back ridge, providing a tight seal and keeping the stopper from backing out. This style is not recommended for vials with a straight-walled finish. The blow-back feature is typically found in vials from North American manufacturers.

TECHNICAL SPECIFICATIONS:

Diameter (mm)	Inner: 7.20 Outer: 13
Thickness (mm)	2.05

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM