



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6149



Image is for illustrative purposes only. Actual product may differ in color, size, and appearance.

FLUID TRANSFER TUBING (STERILE)

Plastic, Male Luer Inlet / Outlet

DESCRIPTION: This is a replacement tubing used in conjunction with the Repeater pump to transfer fluids from one container to another.

FEATURES:

-  Latex free
-  Pyrogen free
-  Sterile when unopened
-  Made from non-DEHP PVC, silicone and ABS