



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6166**



## TUBING ASSEMBLY

For Omnisense and Unispense Peristaltic Pump

**DESCRIPTION:** This tubing assembly is used, in conjunction with the Omnisense and Unispense peristaltic pumps, for various automatic dispensing needs.

### FEATURES:

- Can be sterilized by autoclave or radiation
- Platinum cured silicone translucent tubing, USP Class VI for medical, laboratory and pharmaceutical use

### TECHNICAL SPECIFICATIONS:

Length (m / ft)		15 / 50
Inside Diameter (mm / in)		3.2 / 0.126
Wall Thickness (mm / in)		1.6 / 0.06
Material	Tubing	Silicone
	Sinker	Glass
	Dispensing Nozzle	PP