



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6194



Image is for illustrative purposes only. Actual product may differ in color, size, and appearance.

STIRRING HOT PLATE, CORNING

10" x 10", Digital, 120 V

DESCRIPTION: This equipment allows the user to heat up preparations for melting purposes. In combination with a magnetic stir bar, the preparation can be properly mixed while its ingredients are dissolving. In cases where low heating temperature is required, a water bath can be used to achieve appropriate temperature conditions.

FEATURES:

-  Top plate is easy to clean and highly resistant to scratches, corrosion and chemical attack
-  Twin digital display shows both temperature and stirring speed separately
-  Precise liquid temperatures and stirring speeds
-  Display indicates when set temperature and stirring speeds are reached
-  Safety indicators warns the user when the unit is plugged in and is too hot to touch
-  The angled front panel minimizes the damage from spills
-  Large knobs for better control

WARRANTY: 2 years

TECHNICAL SPECIFICATIONS:

External Dimensions, W x D x H (cm / in)	19.7 x 39.1 x 11.8 / 11.0 x 15.4 x 4.6
Top Plate Dimensions, W x D (cm / in)	25.4 x 25.4 / 10 x 10
Weight (Kg / lb)	5.2 / 11.5
Top Plate Material	Glass-Ceramic
Temperature Range (°C)	5 - 550
Temperature Accuracy (°C)	± 2
Stirring Range (RPM)	60 - 1150
Stirring Speed Accuracy (%)	± 5
Electrical	120 V / 60 Hz / 1113 W / 9.3 A