



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6195



STIRRING HOT PLATE, CORNING

5" x 7", Digital, 120 V

DESCRIPTION: This equipment allows the user to heat up preparations for melting purposes. In combination with a magnetic stir bar, the preparation can be properly mixed while its ingredients are dissolving. In cases where low heating temperature is required, a water bath can be used to achieve appropriate temperature conditions.

FEATURES:

- ✔ Top plate is easy to clean and highly resistant to scratches, corrosion and chemical attack
- ✔ Twin digital display shows both temperature and stirring speed separately
- ✔ Precise liquid temperatures and stirring speeds
- ✔ Display indicates when set temperature and stirring speeds are reached
- ✔ Safety indicators warns the user when the unit is plugged in and is too hot to touch
- ✔ The angled front panel minimizes the damage from spills
- ✔ Large knobs for better control

WARRANTY: 2 years

TECHNICAL SPECIFICATIONS:

External Dimensions, W x D x H (cm / in)	19.7 x 28.0 x 10.8 / 7.8 x 11.0 x 4.3
Top Plate Dimensions, W x D (cm / in)	12.7 x 17.8 / 5.0 x 7.0
Weight (Kg / lb)	3.2 / 7.0
Top Plate Material	Glass-Ceramic
Temperature Range (°C)	5 - 550
Temperature Accuracy (°C)	± 2
Stirring Range (RPM)	60 - 1150
Stirring Speed Accuracy (%)	± 5
Electrical	120 V / 60 Hz / 698 W / 5.9 A