



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6200**




## Vials, Open, Amber Glass, 10 mL, 20 mm (Sterile)

**DESCRIPTION:** All vials are sterile and pyrogen-free. Lot release sterility testing conforms to requirements set forth in USP <71>. Stopper and seal are not included.

### TECHNICAL SPECIFICATIONS:

Filling Volume (mL)	12.5 ± 0.5
Height (mm)	52.0 ± 0.5
Thickness of bottle-lip (mm)	3.9 ± 0.2
Height of bottle-neck (mm)	3.5 ± 0.3
Effective height of bottle-body (mm)	40.0 ± 1.0
Outer diameter of bottle-mouth (mm)	19.8 ± 0.2
Inner diameter of bottle-mouth (mm)	12.6 ± 0.2
Outer diameter of bottle-neck (mm)	15.8 ± 0.3
Body diameter (mm)	22.0 ± 0.2
Material	Borosilicate glass



### APPLICABLE PRODUCTS:

- ☑ 6118: Seals, Aluminum, Center Tear-Out, 13 mm, Silver (Sterile)
- ☑ 6119: Stopper, 20 mm, 2 Prongs, Chlorobutyl, Grey (Sterile)
- ☑ 6141: Seals, Aluminum, Flip-Off, 20 mm, Blank Top, Red (Sterile)
- ☑ 6142: Seals, Aluminum, Flip-Off, 20 mm, Blank Top, Blue (Sterile)
- ☑ 6143: Seals, Aluminum, Flip-Off, 20 mm, Blank Top Green (Sterile)
- ☑ 6145: Stopper, 20 mm, Snap-On, Chlorobutyl, Grey (Sterile)
- ☑ 6146: Stopper, 13 mm, Snap-On, Chlorobutyl, Grey (Sterile)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM