



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6215**







Image is for illustrative purposes only. Actual product may differ in color, size, and appearance.

## **CLEANROOM TACKY MAT, GRAY**

**W 18" × D 46"**

**DESCRIPTION:** This tacky mat helps to remove dirt on contact for cleanroom purposes.

### **FEATURES:**

-  The mat is made up of thirty (30) 2-mil (0.002") thick polyethylene sheets with full width adhesive backing
-  Adhesive is non-allergenic, non-odorous and non-drying and will not age
-  Each sheet contains anti-microbial germicide
-  Each sheet is individually tabbed with a pull-tab system.