



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6230**



## **FREEZONE CONSOLE FREEZE DRYER, LABCONCO (12 L, 230 V)** W 32.4" x D 29.6" x H 63.7", w/Stoppering Tray Dryer

**DESCRIPTION:** This Freeze Dryer is designed for lyophilizing moderate to large sample loads or numerous small container batches. The newly integrated Lyo-Works OS provides real time display of information such as the collector temperature, the vacuum level and the end-point detection system (optional), shown on a 5" full color, capacitive touch screen.

The stoppering tray dryer allows the user to stopper vials or serum bottles under vacuum.

### **FEATURES:**

- ✔ The upright stainless steel collector coil is capable of holding about 12 L of ice before defrosting
- ✔ The console cabinet has a two-piece, removable front panel with viewing window and also an interior space to fit a vacuum pump
- ✔ The chamber lid is made of a 3/4-in thick clear acrylic material and lined with a neoprene gasket
- ✔ The 5" full color capacitive touch screen has the following:
  - ✔ Vacuum set point, alerts, diagnostic tests, user manual, defrost control
  - ✔ Alerts are programmable to send via Ethernet connection to pre-entered email addresses
  - ✔ Maintenance alert settings include vacuum pump oil change, refrigeration cleaning, total operating hours and line voltage offset
  - ✔ Data logging stores and displays data in table or chart format, and may be exported via USB or Ethernet
  - ✔ Level-specific security lock-out prevents changes to programs and other system parameters
- ✔ Vacuum control valve that maintains vacuum level set point
- ✔ Side mounted USB port and rear-mounted Ethernet connection
- ✔ Internal electrical receptacle for vacuum pump connection

**TECHNICAL SPECIFICATIONS:**

Dimensions, W x D x H (cm / in)	82.2 x 75.1 x 161.7 / 32.4 x 29.6 x 63.7
Shipping Weight (kg / lb)	337 / 743
Electrical (Console and Tray Dryer)	230 V, 60 Hz, 9 A
Plug Type	North America plug
Vacuum Pump Displacement Required	163 L/min or larger
Collector Temperature (°C / °F)	-84 / -119
Ice Holding Capacity (L)	12
Conformance	CAN/CSA C22.2, UL

**WARRANTY:** 1 year