



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6279






AIR FLOW TEST TUBE KIT


DESCRIPTION: This air flow test tube kit can be used universally. It is simple and easy to use for the identification of air flow.

Included in this kit: Plastic carrying case (1), Tube Breaker (1), Air Flow Test Tubes (10), Rubber Caps (10), Rubber Inflation Bulb, and related documentation.

FEATURES:

-  Quick information: Provides source, direction and speed of air flow immediately
-  All-purpose: The test tube can be used until the smoke has completely disappeared. If the test is completed before all the smoke disappears, the tube can be sealed again using the included rubber cap
-  Simple to use: Smoking sulfuric acid is located in the tube. When the top of the tube is opened, a small rubber ball is used to pump air through the tube. This creates a visible white smoke, which is carried on any existing flow.

ADDITIONAL INFORMATION:

-  Product No. 6461 are additional test tubes which can be used with the accessories of this kit.