




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6299



DISPERSING ELEMENT, IKA (0.5 – 20 L) For T 50 Digital Dispersing Instrument

DESCRIPTION: This dispersing element is an accessory to the T 50 Dispersing instrument. It functions to properly homogenize and emulsify solutions and prepare good suspensions.

TECHNICAL SPECIFICATIONS:

Diameter, Stator / Rotor (mm)	45 / 36
Shaft Length (mm)	300
Immersion Depth (mm)	70 – 250
Working Range, Water (L)	0.5 – 20
Material of Construction	PTFE and AISI 316L
Maximum Allowable Speed (RPM)	10000
Maximum Temperature (°C)	180
Ultimate Fineness, Suspensions (µm)	40 – 100
Ultimate Fineness, Emulsions (µm)	10 – 30

WARRANTY: 2 years