



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.:



SAMIX MIXING TOOLS

Disposable Mixing Blades and Shafts

DESCRIPTION: The disposable mixing blades provide a trituration effect and improved mixing. The disc combines with the reusable mixing shafts to form the disposable mixing tool. The disc can be removed after each preparation or can be left inside the jar.

AVAILABLE SIZES:

See the tables below for the available sizes for the disposable blades and the **required** shafts:

Table 1: Disposable blades

Capacity (mL)	Product No.
15 – 30	6303
50	6302
100 – 200	6301

Table 2: Required blade shafts

Size	Product No.
15 – 100 ml	6649
200 ml	6648
Set (15 – 100 mL [1], 200 mL [1])	6647