



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6387



## RPCR (REMOTE PC REPORTING) SOFTWARE, TUTTNAUER

For EZPlus Autoclaves

**DESCRIPTION:** This software generates a report that includes detailed specifications of the cycle including graphs and tables to measure specific parameters and data values, with an option to save the report in PDF file format. This software can be operated either in online mode, via Ethernet connection, or offline mode via USB.

### FEATURES:

- Manage the history report files of the cycles run on the autoclave
- Track all the history of the processes that have been run on the autoclave
- Track which parameters have been used in each of the cycles and stages used
- Choose the type of the report between graph, table or print out
- Generate a report in PDF format, and choose the type of report
- The graph type application can offer the user an option to add some more inputs or outputs to the presentation