



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6392



PROFILLER 3700 CAPSULE FILLING SYSTEM, TORPAC (Size 5) 300-Hole Capsule Machine

DESCRIPTION: The ProFiller 3700 is the new standard in 300-hole bench-top capsule filling systems. Experienced users of the profiling system benefit greatly from the higher output, fast size changes, easy to clean design and optional accessories. First time users like that it is easy to use and the different upgrade options.

Included Parts: Hand-Held Filler Base, Orienter Base, Locking Plate, 60-Pin Tamper, Powder Spreader, Powder Tray, Powder Tray Clamps, Plastic Storage Box, Capsule Bins, Measuring Bowls and Change Parts (Caps Tray, Body Sheet Set (A, B, C), Spacer Plate, Capsule Locker and Orienter Tray).

FEATURES:

- ✔ Shipped fully assembled. Ready for use in less than 1 hour
- ✔ Easier to use than other 200 or 300-hole systems
- ✔ All parts are dishwasher safe. Can be dried in a hot air oven or with a portable hair dryer
- ✔ Change capsule sizes in only 5 minutes
- ✔ Orienter Base accommodates capsule sizes 00 – 4
- ✔ All metal parts Inox (SS316) for reduced static
- ✔ Unique features include lightweight 2-piece Orienter, Orienter Locator and Oversize Powder Tray & Push Plate
- ✔ Designed for fast disassembly & visual verification of cleaning
- ✔ Fill weights, USP/EP method within $\pm 1-3\%$ of average
- ✔ Works with all capsules including gelatin, vegetarian/HPMC and starch capsules
- ✔ Tamper size range optimized for higher fill weights and faster filling
- ✔ Indicator gauges to assist in capsule orientation, separation and locking
- ✔ Fits under most powder hoods and in glove boxes
- ✔ 3700 fully locks capsules after filling, unlike most hand-held fillers
- ✔ Compliance: CE, cGMP

WARRANTY: 1 year