



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6401



WASTE CHUTE BAG, LABCONCO

Clear, For XPert and RXPert Enclosures

DESCRIPTION: This bag attaches to the waste port and allows for temporarily holding waste while working under the enclosure.

TECHNICAL SPECIFICATIONS:

Dimensions, W × L (cm / in)	25.4 × 40.6 / 10 × 16
Shipping Weight (kg / lb)	2.3 / 5
Thickness (in)	0.004