



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6406









CIMAREC+ STIRRING HOT PLATE, THERMO SCIENTIFIC

7.25" x 7.25", Digital

DESCRIPTION: This digital stirring hotplate offers precise stirring control, safety and temperature performance for routine protocols.

FEATURES:

-  Slow-speed stirring, consistent speed control, and strong magnetic coupling
-  HOT TOP warning system to protect from accidental burns with prominent display when heating surface is above 50 °C (122 °F)
-  Automatic over-temperature protection
-  Easy-to-read LED display for heating and stirring
-  Raised display design protects electronics from spills
-  Stir bar included

WARRANTY: 3 years

TECHNICAL SPECIFICATIONS:

External Dimensions, W x D x H (cm / in)	20.8 x 36.0 x 10.5 / 8.2 x 14.2 x 4.1
Top Plate Dimensions, W x D (cm / in)	18.4 x 18.4 / 7.25 x 7.25
Weight (Kg / lb)	3.4 / 7.5
Top Plate Material	Aluminum
Temperature Range (°C / °F)	30 – 300 (86 – 572)
Stirring Range (RPM)	50 – 1500
Conformity	EN 61000-3-2 EN 61000-3-3 EN 61326-1 EN 61010-1 EN 61010-2-010 EN 61010-2-051 CE (230 V Models)