



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6409



HOSE BARB TO FEMALE LUER CONNECTOR

Stainless Steel, 2 – 3 mm (Inside Diameter) Tubes

DESCRIPTION: This luer adapter allows connection to equipment and components with a female luer connector.

TECHNICAL SPECIFICATIONS:

Inside Diameter (mm / in)	2 / 0.07
Length (mm / in)	3 / 0.108
Material	Stainless Steel
For Use With	Dispensing Pump
Type	Hose Barb/Luer Adapter