



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6410**



**PLASTIC OINTMENT TUBE w/SCREW CAP (8 oz)**  
4.7 mm Tube Opening, White

**DESCRIPTION:** This plastic ointment tube is lightweight yet durable to withstand squeezing and folding without developing cracks or leaks. The tube is made from FDA compliant Low Density PolyEthylene (LDPE) and should be heat sealed to close.

**TECHNICAL SPECIFICATIONS:**

Dimensions w/Cap, Ø x L (in)	2 x 7
Capacity (mL / oz)	250 / 8
Material	LDPE
Exterior color	White
Liner	None
Orifice (mm / in)	4.7 / 0.187