



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6411**







## SELF SEAL STERILIZATION POUCH (9.5" × 13")

Nylon, Dry Heat

**DESCRIPTION:** This pouch is manufactured for use with dry heat processes up to 420 °F. The nylon material is clear, strong, and puncture resistant, maintaining a durable barrier for instrumentation until point of use.

### FEATURES:

-  Latex-free
-  Built-in adhesive tape which eliminates the need for a heat sealer or indicator tape
-  Pouch notched for ease of opening after processing
-  Complimentary external indicator labels included

### TECHNICAL SPECIFICATIONS:

Dimensions (in)	9.5 × 13
Material	Nylon
Max. Temperature (°F)	420