



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6419



REFRIGERATOR (Left Hand Door Swing, 481 L / 17 cu. ft)
W 28" x D 25" x H 80", 36 °F – 47 °F / 2 °C – 8 °C, Auto Defrost

DESCRIPTION: This is a full-size all-refrigerator designed for optimum pharmaceutical storage of sensitive materials. This environmentally friendly refrigerator comes with a glass door and stainless steel interior and exterior finish.

FEATURES:

- ✔ Commercially approved: ETL-S listed to ANSI-NSF Standard 7 and meets UL-471
- ✔ Large capacity: 17 cu. ft interior for maximum storage and convenience
- ✔ LED lighting: Individual lights in door pane offer a full view of stored contents, with an off switch located inside
- ✔ Automatic defrost: Reduced user maintenance with automatic defrost
- ✔ Stainless steel interior and exterior finish: Cabinet and interior walls are constructed from stainless steel
- ✔ Glass door: Tempered glass door provides heat-safe view of interior contents
- ✔ Factory installed lock: Keyed lock for a secure interior
- ✔ Adjustable cantilevered shelves: 6 molded plastic shelves can be arranged in virtually any position to accommodate a variety of sizes
- ✔ Self-closing door: Door automatically closes if left ajar to ensure items are not accidentally overexposed to warm temperatures

TECHNICAL SPECIFICATIONS:

Dimensions, W x D x H (cm / in)	70 x 63 x 203 / 27.63 x 24.75 x 79.75
Depth with door at 90° (cm / in)	130 / 51
Capacity (L / cu. ft)	481 / 17
Weight	95 / 210
Shipping Weight (kg / lb)	118 / 260
Electrical	115 V AC / 60 Hz / 5 A
Cabinet	Stainless Steel
Shelf type	Plastic
Shelf Quantity	6
Adjustable Shelves	Yes
Thermostat Type	Digital
High Side Pressure (PSI)	186
Low Side Pressure (PSI)	86

Interior Light	Yes
Temperature Range (°C / °F)	2 – 8 / 36 – 47
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Sanitation	ETL-S

WARRANTY: 1 year