



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6496**



**SHARPSAFETY RCRA HAZARDOUS WASTE CONTAINER, COVIDIEN (18 Gallon)**  
Polypropylene, Black w/Hinged lid

**DESCRIPTION:** The RCRA Hazardous Waste Containers are designed to segregate hazardous pharmaceutical waste under the Resource Conservation and Recovery Act (RCRA). Black with lids, they are easily distinguished from other disposal containers.

**FEATURES:**

- ✔ RCRA Hazardous waste label meeting EPA satellite accumulation labeling requirements
- ✔ White and black color makes it easily distinguished from sharps, chemo and other disposal containers
- ✔ Toxic and flammable additional labels included to designate contents
- ✔ Tamper evident label helps deter tampering during final disposal
- ✔ PGII Rated. Meets DOT requirements for transport without secondary packaging
- ✔ Leak-resistant gasket and absorbent pad helps to contain liquid contents should a primary container break
- ✔ The containers can be nested to save valuable storage space