



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6578







3D INVERSION MIXING ARM

For Lab 1 Universal Mixer

DESCRIPTION: This is an attachment for the Lab 1 Universal Mixer. This arm should be used when in-container mixing of liquids or powders is required.

FEATURES:

-  3D twin inversion motion for quick efficient mixing
-  In container mixing
-  Can be used to mix liquids or powder
-  Fits vessels up to 750 mL