



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6579










UNIVERSAL LAB MIXER – THE LAB 1, BASE UNIT ONLY

W 9" x D 11" x H 8"

DESCRIPTION: The Lab-1 universal mixer allows the user to obtain a perfect homogenized mix for both liquids and powder. The Lab-1 is also designed for “in-container” mixing of any size and shape.




FEATURES:

-  Unique twin inversion motion for quick efficient mixing
-  In container mixing
-  Small, quiet and compact
-  Speed control from 0 – 60 rpm
-  Interchangeable mixing arms for greater flexibility
-  Adjustable angle from horizontal through 45 degrees
-  Upgrade path for future attachments

TECHNICAL SPECIFICATION:

Dimensions (base unit only), W x D x H (mm / in)	229 x 290 x 203 / 9 x 11 x 8
Weight, kg / lb (base unit only)	8.0 / 17.6
Speed Range (rpm)	0 – 60
Electrical Requirements	110 – 240 V, 50 / 60 Hz
Maximum container size, D x H (mm / in)	101 x 127 / 4 x 5

RELATED ACCESSORIES

-  6578-01: 3-dimensional inversion mixing arm
-  6577-01: V-shell mixing container
-  6576-01: Tube rotation disc, for multiple vial mixing

WARRANTY: 1 year