



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6581



BOTTLE HOLDER, LABCONCO
Steel, epoxy-coated

DESCRIPTION: This epoxy-coated steel bottle holder attaches to the rear baffle perforations to hold a squirt bottle or other round container. It also helps to maximize the usable work surface within XPert Enclosures, Stations and Systems to make the workspace more efficient.

TECHNICAL SPECIFICATION:

Opening, Ø (mm / in)	76.2 / 3.0
Weight (kg / lb)	0.5 / 1.0