



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6595



Mixer, Samix U1000

DESCRIPTION: The Samix machine is used to efficiently mix creams, ointments and semi-solids. It allows to incorporate and mix all necessary ingredients, in their specific measurements, in the jar saving lots of processing time. The U1000 can be used to mix preparations in the 15 – 1000 mL jars.

FEATURES:

- ✔ Automatic or manual lifting function
- ✔ 10-steps (electronically controlled) stirring speed regulation
- ✔ Variable time setting, from 5 seconds to 9:55 (9 minutes and 55 seconds)
- ✔ For preparations up to 1000 mL
- ✔ Includes 15-30mL, 50mL, 100mL, 200mL, 300mL, 500mL, and 1000mL standard mixing blades
- ✔ Power cord included

TECHNICAL SPECIFICATIONS:

Dimension, W × D × H (cm / in)	22.1 × 35.6 × 64.2 / 8.70 × 14.0 × 25.3
Weight (kg / lb)	16.7 / 36.8
Stirrer Motor Power (W)	550
Lifting Motor Power (W)	50
Electrical Requirements	100 – 240 V / 50 – 60 Hz
Speed (RPM)	250 – 2500
Protection Type	IP X0
PC connection	USB - Standard B
Protection Class	I
Operation Class	S1
Conformity	CE

WARRANTY: 2 years

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM