



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6597**



## Mixer, Samix ES500

**DESCRIPTION:** The Samix machine is used to efficiently mix creams, ointments and semi-solids. It allows to incorporate and mix all necessary ingredients, in their specific measurements, in the jar saving lots of processing time. The ES500 can be used to mix preparations in the 15 – 500 mL jars.

### FEATURES:

- ✔ Automatic or manual lifting function
- ✔ 10-steps (electronically controlled) stirring speed regulation
- ✔ Variable time setting, from 5 seconds to 9:55 (9 minutes and 55 seconds)
- ✔ Maximum operation time is 5 minutes under high load
- ✔ For preparations up to 500 mL
- ✔ Includes 15-30mL, 50mL, 100mL, 200mL, 300mL, and 500mL standard mixing blades
- ✔ Power cord included

### TECHNICAL SPECIFICATIONS:

Dimension, W × D × H (cm / in)	30 × 30 × 65 / 11.8 × 11.8 × 25.6
Weight (kg / lb)	9.5 / 20.9
Stirrer Motor Power (W)	220
Lifting Motor Power (W)	50
Electrical Requirements	230 V 50 Hz, 100/120 V 60 Hz
Speed (RPM)	650 – 2100
Protection Type	IP 21
Protection Class	I
Operation Class	KB5
Conformity	CE

**WARRANTY:** 2 years

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)