



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING


## PRODUCT NO.:




## SAMIX SPINDLE

**DESCRIPTION:** The spindle is an extracting tool that is used in combination with some jars to make dispensing easy. The spindle also comes installed in some White Jar/Lid combinations, see Table 2 below.

*Table 1: Spindle only*

|               | <br>Spindle Only |
|---------------|---|
| Jar Size (mL) | Product No.   |
| 300           | 6603  |
| 500           | 6602  |

*Table 2: Jar w/Spindle and Lid*

|               |  |
|---------------|--|
| Jar Size (mL) | Product No.  |
| 300           | 6628   |
| 500           | 6626   |